



## Security with dignity. Safety with discretion.

*Wander management solutions from Ascom teleCARE IP*



- Discreet monitoring and alert solution that promotes residents' dignity, autonomy and self-esteem
- Enables mobility, exercise and socialization in a safe, controlled environment
- Reassures residents' families; provides data and records for audits; helps ensure care quality
- Accurate location of wandering/residents supports swift staff responses. Loitering alerts are triggered when at-risk residents are in the vicinity of certain exits beyond a preset time.
- Helps minimize intrusive monitoring. Frees time for more meaningful staff-resident interactions

Ascom teleCARE IP is a proven way to reliably and discreetly manage wander and elopement behavior in long-term elderly care facilities. In use at thousands of facilities worldwide, teleCARE IP deploys wearable transceivers and low-frequency (LF) beacons to automatically alert staff or lock exits in cases of elopement and loitering. Such functionality is only part of a range of options that can be fully integrated into a complete emergency and nurse call system.

The transceiver worn by residents can either be worn as a pendant, or as a wristband that can be fastened when needed. The latter features a tamper alarm, and is often preferred for those with advanced dementia. Both versions can be color-coded for easy resident recognition, and can include a 'no-motion' feature that alerts staff when no movement has been detected over a pre-determined period of time.

### **Part of an integrated solution**

Wander, elopement and loiter management functionality is an integral part of an Ascom teleCARE IP nurse call solution. It is not a separate component or installation. It does not require separate management systems or dedicated mobile devices. Wander alerts can be sent to the same device—an Ascom Myco 3 smartphone, for instance—that staff use for other communication, coordination and collaboration tasks.

Transmitting wander alerts to caregivers' mobile devices has the added benefit of helping to foster a calmer, quieter environment. This is obviously beneficial for staff at risk of alarm fatigue. It is also crucial for residents' well-being, as loud and disruptive electronic alerts can cause the dementia resident to become anxious or agitated—mental states known to trigger wandering behavior.<sup>5</sup>



## Customized and automated alerts and actions


An Ascom teleCARE IP solution is highly customizable. Functionalities and features can be tailored to match existing and/or planned infrastructures and systems. Wander and loitering management functions, for example, can be configured to match any type/size/layout of aged-care facility—and of course the specific needs of individual residents.

### *Dementia-related wandering and elopement*

**Wandering** is defined by the North American Nursing Diagnosis Association as: “meandering, aimless, or repetitive locomotion that exposes the individual to harm; frequently incongruent with boundaries, limits or obstacles.”<sup>1</sup>

**Elopement** can be understood as “the act of leaving a safe area unsupervised and unnoticed and entering into harm’s way.”<sup>2</sup> According to the Alzheimer’s Association, six out of ten dementia sufferers will wander, putting themselves at risk of injury and/or becoming lost.<sup>3</sup>

It is estimated that almost **47 million people** currently live with dementia, with that number rising to 75 million by 2030.<sup>4</sup>



Wander alerts can be sent to the same device—an Ascom Myco 3 smartphone, for instance—that staff use for other communication, coordination and collaboration tasks. Transmitting wander alerts to caregivers’ mobile devices has the added benefit of helping to foster a calmer, quieter environment.

For Sales and Support please contact your local Ascom partner



Connexus Communications Limited  
Unit D2  
Baldonnell Business Park  
Naas Road  
Dublin 22  
Ireland  
T: 00 353 1 4105866, F: 00 353 1 4105867  
info@connexus.ie  
www.connexus.ie

#### Ascom Holding AG

Zugerstrasse 32,  
CH-6340 Baar,  
Switzerland  
info@ascom.com  
Phone: +41 41 544 78 00  
ascom.com

**ascom**

1. Miller-Keane, O’Toole, M.T. (2006) Miller-Kean Encyclopedia and Dictionary of Medicine, Nursing, and Allied Health, 7th edn. Philadelphia: Saunders <<http://medical-dictionary.thefreedictionary.com/wandering>> [20 May 2019]
2. Moishita, L. (1990) ‘Wandering Behavior’, in Alzheimer’s Disease: Treatment and Long-Term Management, eds. Cummings, J.L., Miller B.L. New York: Dekker, 157-176
3. Alzheimer’s Association (2016) Wandering and Getting Lost [online] available from <<http://www.alz.org/care/alzheimers-dementia-wandering.asp>> [20 May 2019]
4. Prince, M., Wimo, A., Guerchet, M., Ali, G., Wu, Y., Prina, M. (2015) World Alzheimer Report 2015: The Global Impact Of Dementia [online] London: Alzheimer’s Disease International. Available from <<http://www.alz.co.uk/research/WorldAlzheimerReport2015.pdf>> [20 May 2019]
5. Alzheimer’s Australia (2016) Wandering [online] available from <<https://fightdementia.org.au/national/support-and-services/carers/behaviour-changes/>> [20 May 2019]

Multiple products from Ascom and third parties support the features and functions described. All Ascom solutions may not be available in specific countries.